

# ODM98-MD

<http://www.gigahertz-optik.de/de-de/produkt/ODM98-MD>

**Produkt-Tags:**



# Überblick

ODM98 is commonly used as reflectance standards (flat surface) for the calibration of:

- Spectrophotometers
- Imaging System Pixel Uniformity
- Reflectance Meters
- Surface Color Meters

**ODM98-MD** are plates in disc shape with surfaces machined to a flat optical finish. Gigahertz-Optik GmbH stocks different sizes of 10.5mm thick discs. For OEM applications custom sizes, shapes and thicknesses are supplied.

**BN-R98** Reflectance Standards are ODM98-MP plates or ODM98-MD discs supplied in a protective housing with removable cover. Gigahertz-Optik GmbH stocks different sizes of BN-R98 reflectance standards in disc and plate shape. All standards can be supplied with optional spectral reflectance calibration from 250 nm to 2400 nm in 50nm steps. Calibration is done with a double-beam spectrophotometer in 8/d measurement geometry using reflectance standards traceable to national calibration laboratories. For OEM applications custom sizes, shapes and housings are available.

## Technische Daten

Spezifikationen	
ODM98-B rohes Plattenmaterial	NULL
ODM98-B rohes Folienmaterial	NULL
ODM98-B rohes Blockmaterial	Abmessungen
ODM98-B01	200 x 200 x 200 mm
ODM98-B02	120 x 120 x 250 mm
ODM98-B03	120 x 120 x 120 mm
ODM98-B04	65 x 65 x 500 mm
ODM98-B05	65 x 65 x 245 mm
ODM98-B06	65 x 65 x 120 mm
ODM98-B07	D170 x 190mm
ODM98-B08	76 x 76 x 250 mm
ODM98-B13	D52 x 250 mm
ODM98-P01	560 x 260 x 12 mm
ODM98-P02	305 x 305 x 12 mm

ODM98-P03	260 x 260 x12 mm
ODM98-P04	135 x 135 x 12 mm
ODM98-P05	55 x 55 x 12 mm
ODM98-F01	1.5 x 300 x 1000 mm*
ODM98-F02	0.2 x 105 x 1000 mm*
ODM98-F03	0.5 x 105 x 1000 mm*
NULL	*) minimum length cut from roll

## Bestellinformationen

Artikel-Nr	Modell	Beschreibung
<b>Produkt</b>		
15295954	ODM98-B01	ODM98-B Raw Block Material
15295955	ODM98-B02	ODM98-B Raw Block Material
15295956	ODM98-B03	ODM98-B Raw Block Material
15295957	ODM98-B04	ODM98-B Raw Block Material
15295958	ODM98-B05	ODM98-B Raw Block Material
15295959	ODM98-B06	ODM98-B Raw Block Material
15295960	ODM98-B07	ODM98-B Raw Block Material
15296490	ODM98-B08	ODM98-B Raw Block Material
15296341	ODM98-B13	ODM98-B Raw Block Material
15296491	ODM98-P01	ODM98-B Raw Plate Material
15296492	ODM98-P02	ODM98-B Raw Plate Material
15296493	ODM98-P03	ODM98-B Raw Plate Material
15296494	ODM98-P04	ODM98-B Raw Plate Material
15296495	ODM98-P05	ODM98-B Raw Plate Material
15296022	ODM98-F01	ODM98-B Raw Foil Material
15295963	ODM98-F02	ODM98-B Raw Foil Material
15296241	ODM98-F03	ODM98-B Raw Foil Material